

World Standard for blood pressure

The Disposable Calibrated® V-Lok® Cuff

Our **Disposable Calibrated® V-Lok® Cuffs** are *latex-free* and are fully interchangeable with our complete line of reusable cuffs. They are intended for **single patient use** to help prevent cross contamination. All sizes are available with one or two tubes and a variety of connectors for use with manual sphygmomanometers or electronic non-invasive blood pressure monitors.

As with all Baumanometer® Calibrated® V-Lok® Cuffs our **Disposable Calibrated® V-Lok® Cuffs** feature:

- The Calibrated® V-Lok® system which immediately and clearly indicates whether or not the correct size is being used.
- Genuine Velcro® hook and loop fasteners which are precisely positioned and evenly matched in holding strength to provide the optimal number of open-close cycles.
- Cuffs that clearly indicate the center of the bladder on the inner and outer surfaces for ease of application on either right or left limb.
- The high quality design, materials, and workmanship of products made in the USA.

The Disposable Calibrated® V-Lok® Cuff - a World Standard . . . in every detail.


Available through better-known medical equipment dealers worldwide.



W.A. Baum Co., Inc.

620 Oak Street
Copiague, NY 11726 USA
www.wabaum.com

Tel: 631-226-3940
Fax: 631-226-3969
info@wabaum.com



Disposable Calibrated® V-Lok® Cuffs

Ordering Information:

Disposable Calibrated® V-Lok® Cuffs - Double Tube (1 Female Luer Connector)

These cuffs are a latex-free alternative for our reusable cuffs.

Cat. No.	Description	Limb Circumference
5260	Adult	25 cm to 35 cm
5261	Child/Small Adult	18 cm to 26 cm
5262	Infant	10 cm to 19 cm
5265	Large Arm	33 cm to 47 cm

Disposable Calibrated® V-Lok® Cuffs - Double Tube (2 Threaded Adapters)

Each cuff has two threaded adapters for use with Dinamap® monitors.

Cat. No.	Description	Limb Circumference
5260TA	Adult	25 cm to 35 cm
5261TA	Child/Small Adult	18 cm to 26 cm
5262TA	Infant	10 cm to 19 cm
5265TA	Large Arm	33 cm to 47 cm

Disposable Calibrated® V-Lok® Cuffs - Single Tube (No Connector)

Cat. No.	Description	Limb Circumference
5130	Adult	25 cm to 35 cm
5131	Child/Small Adult	18 cm to 26 cm
5132	Infant	10 cm to 19 cm
5135	Large Arm	33 cm to 47 cm

Disposable Calibrated® V-Lok® Cuffs - Single Tube (Female Luer Lock)

Cat. No.	Description	Limb Circumference
5130LL	Adult	25 cm to 35 cm
5131LL	Child/Small Adult	18 cm to 26 cm
5132LL	Infant	10 cm to 19 cm
5135LL	Large Arm	33 cm to 47 cm

Disposable Calibrated® V-Lok® Cuffs - Single Tube (CPC Connector)

These cuffs are used with PROTOCOL® PROPAQ monitors.

Cat. No.	Description	Limb Circumference
5130CPC	Adult	25 cm to 35 cm
5131CPC	Child/Small Adult	18 cm to 26 cm
5132CPC	Infant	10 cm to 19 cm
5135CPC	Large Arm	33 cm to 47 cm

Disposable Calibrated® V-Lok® Cuffs - Single Tube (Plastic Male Rectus Connector)

These cuffs are used with SIEMENS MEDICAL SYSTEMS and other non-invasive blood pressure monitors.

Cat. No.	Description	Limb Circumference
5130OC	Adult	25 cm to 35 cm
5131OC	Child/Small Adult	18 cm to 26 cm
5132OC	Infant	10 cm to 19 cm
5135OC	Large Arm	33 cm to 47 cm

Disposable Calibrated® V-Lok® Cuffs - Single Tube (1M/2F Luer Locks)

These cuffs are used with SPACE LABS monitors. They are equipped with one male luer lock, two female luer locks and nine feet of straight tubing.

Cat. No.	Description	Limb Circumference
5130SPACE	Adult	25 cm to 35 cm
5131SPACE	Child/Small Adult	18 cm to 26 cm
5132SPACE	Infant	10 cm to 19 cm
5135SPACE	Large Arm	33 cm to 47 cm

Catalog Number Guide: TA = Threaded Adapter, CPC = Colder Products Connector, OC = Rectus Connector, LL = Luer Lock which are either M = Male or F = Female.

Note: Baumanometer, Standby, V-Lok, Extendex, Calibrated, Hi/Lo, and Air-Flo are registered trademarks of W.A. Baum Co., Inc.

W.A. Baum Co., Inc.

620 Oak Street
Copiapue, NY 11726 USA
www.wabaum.com

Tel: 631-226-3940
Fax: 631-226-3969
info@wabaum.com

Since 1916 Originator and Maker of Blood Pressure Apparatus Exclusively

Baumanometer®